


# pDHB1-largeT Control Vector P03503



## Vector features

Position	Feature	Position	Feature
Start: 107 End: 1580	<i>ADH1</i> promoter	Start: 4632 End: 4888	<i>CYC1</i> terminator
Start: 1598 End: 1705	<i>Ost4</i>	Start: 5621 End: 7831	<i>LEU2</i> auxotrophic marker
Start: 1721 End: 3595	<i>largeT</i> antigen fragment	Start: 7980 End: 9013	<i>KanR</i> resistance gene
Start: 3611 End: 3745	<i>Cub</i>	Start: 9306 End: 9973	<i>pBS</i> replication origin ( <i>E. coli</i> )
Start: 3770 End: 4375	<i>LexA</i>	Start: 10095 End: 10634	<i>CEN/ARS</i> replication origin (yeast)
Start: 4393 End: 4630	<i>VP16</i>	Start: 10747 End: 51	<i>lacP</i>