


pDSL-Δp53 Control Vector P03504



Vector features

Position	Feature	Position	Feature
Start: 112 End: 402	CYC1 promoter	Start: 2568 End: 3761	TRP1 auxotrophic marker
Start: 414 End: 530	NubG, amino acids 1-38 of ubiquitin	Start: 3825 End: 5172	2micron replication origin (yeast)
Start: 531 End: 560	HA epitope tag	Start: 6298 End: 6965	pBS replication origin (<i>E. coli</i>)
Start: 603 End: 1538	Delta-p53, N-terminally truncated human p53	Start: 5305 End: 6163	Amp, ampicillin resistance gene
Start: 1592 End: 1853	CYC1 terminator	Start: 7185 End: 50	lacP