


pACT-largeT control vector P03114



Vector features

Position	Feature	Position	Feature
Start: 7392 End: 7789	ADH1 promoter	Start: 4415 End: 4742	ADH1 terminator
Start: 7314 End: 7376	SV40 nuclear localization signal	Start: 2474 End: 3568	LEU2 open reading frame
Start: 6957 End: 7313	GAL4 activation domain	Start: 1 End: 2055	2micron origin of replication
Start: 6929 End: 6956	HA epitope tag	Start: 9015 End: 9945	AmpR resistance gene
Start: 5019 End: 6887	SV40 largeT antigen, amino acids 87-708	Start: 8224 End: 8867	pBR322 origin of replication