



P03103

pGAD-HA library vector



Multiple cloning site



Vector features

Position	Feature	Position	Feature
Start: 22 End: 1476	ADH1 promoter	Start: 2739 End: 3884	LEU2 auxotrophic marker
Start: 1522 End: 1542	SV40 nuclear localization signal	Start: 4651 End: 5487	pUC origin of replication
Start: 1492 End: 1899	GAL4 activation domain	Start: 5644 End: 6504	AmpR resistance gene
Start: 1936 End: 1962	HA epitope tag	Start: 7068 End: 8053	2micron origin of replication
Start: 2349 End: 2675	ADH1 terminator		