


pLexA-p53 control vector P03112



Vector features

Position	Feature	Position	Feature
Start: 71 End: 414	<i>ADH1</i> promoter	Start: 4214 End: 5204	KanR resistance gene
Start: 433 End: 1038	LexA	Start: 5226 End: 6593	2micron origin of replication
Start: 1060 End: 2025	human p53, amino acids 73-393		
Start: 2271 End: 3096	<i>TRP1</i> auxotrophic marker		
Start: 3266 End: 4031	<i>colE1</i> origin of replication		