


pDL2-Alg5 Control Vector



Vector features

Position	Feature	Position	Feature
Start: 112 End: 821	<i>ADH1</i> promoter	Start: 3008 End: 4201	<i>TRP1</i> auxotrophic marker
Start: 842 End: 1843	Alg5 open reading frame	Start: 4265 End: 5612	2micron origin of replication
Start: 1868 End: 1900	HA epitope tag	Start: 5745 End: 6603	AmpR resistance gene
Start: 1907 End: 2044	NubG, amino acids 1-38 of ubiquitin	Start: 6738 End: 7405	pBluescript origin of replication
Start: 2037 End: 2294	<i>CYC1</i> terminator		