


P03242

pNubG-Fe65 prey construct



Position	Feature
Start: 343 End: 1064	ADH1 promoter
Start: 1065 End: 1178	NubG, ubiquitin amino acids 1-38
Start: 1188 End: 1221	HA epitope tag
Start: 1238 End: 3237	Rat Fe65 coding sequence
Start: 3273 End: 3469	ADH1 terminator
Start: 3828 End: 4495	pBluescript origin of replication
Start: 4646 End: 5530	Ampicillin resistance gene
Start: 5621 End: 6968	2micron origin of replication
Start: 7517 End: 8188	TRP1 auxotrophic marker

