


P03240

pOst1-NubI control vector



Position	Feature
Start: 60 End: 767	<i>ADH1</i> promoter
Start: 770 End: 2302	<i>Saccharomyces cerevisiae</i> Ost1 coding sequence
Start: 2303 End: 2414	NubI, amino acids 1-38
Start: 2460 End: 2652	<i>ADH1</i> terminator
Start: 2660 End: 3831	pUC origin of replication
Start: 4691 End: 3834	Ampicillin resistance gene
Start: 4822 End: 6169	2micron origin of replication
Start: 6705 End: 7376	<i>TRP1</i> auxotrophic marker