



pPR3-STE DUALmembrane System 3 Prey Vector

P03237



Vector features

Position	Feature	Position	Feature
Start: 63 End: 770	ADH1 promoter	Start: 2048 End: 3241	TRP1 auxotrophic marker
Start: 788 End: 835	STE2 leader sequence	Start: 3305 End: 4652	2micron origin of replication
Start: 908 End: 940	HA epitope tag	Start: 4785 End: 5643	AmpR resistance gene
Start: 947 End: 1081	NubG, ubiquitin amino acids 1-38	Start: 5778 End: 6445	pBluescript origin of replication
Start: 1077 End: 1334	CYC1 terminator		