


pKS1-ST-Gal control vector

Vector features

Position	Feature	Position	Feature
Start: 68 End: 630	<i>ADH2</i> promoter	Start: 5184 End: 6187	KanMX cassette
Start: 641 End: 667	Strep tag	Start: 6384 End: 7751	2micron replication origin (yeast)
Start: 674 End: 703	HA tag	Start: 8107 End: 8967	AmpR resistance gene
Start: 710 End: 3952	<i>lacZ</i> (<i>E. coli</i> β -galactosidase)	Start: 9115 End: 9782	pBS replication origin (<i>E. coli</i>)
Start: 3954 End: 4215	<i>CYC1</i> terminator		