


pKS1-ST Strep-tag expression vector P03306



```

          BamH I   Sma I               EcoR I               Sal I   Xho I
          -----
cca gat tac gct gga tcc ccc ggg ctg cag gaa ttc gat atc aag ctt atc gat acc gtc gac ctc gag
P   D   Y   A   G   S   P   G   L   Q   E   F   D   I   K   L   I   D   T   V   D   L   E
----->
Strep/HA-tag

```

Vector features

Position	Feature	Position	Feature
Start: 68 End: 630	ADH2 promoter	Start: 1986 End: 2989	KanMX cassette
Start: 641 End: 667	Strep-tag	Start: 3186 End: 4553	2micron replication origin (yeast)
Start: 674 End: 703	HA tag	Start: 4909 End: 5769	AmpR resistance gene
Start: 704 End: 761	Multiple cloning site	Start: 5917 End: 6584	pBS replication origin (<i>E. coli</i>)
Start: 756 End: 1017	CYC1 terminator		