


pKS2-ST Strep-tag secretion vector P03308



```

BamH I   Sma I           EcoR I           Sal I   Xho I
-----
gct gga tcc ccc ggg ctg cag gaa ttc gat atc aag ctt atc gat acc gtc gac ctc gag
 A  G  S  P  G  L  Q   E  F  D  I  K  L  I  D  T  V  D  L  E
----->
SUC2-Strep-HA reading frame

```

Vector features

Position	Feature	Position	Feature
Start: 68 End: 630	ADH2 promoter	Start: 822 End: 1083	CYC1 terminator
Start: 641 End: 706	SUC2 export signal	Start: 2052 End: 3055	KanMX cassette
Start: 707 End: 733	Strep-tag	Start: 3252 End: 4619	2micron replication origin (yeast)
Start: 740 End: 769	HA tag	Start: 4975 End: 5835	AmpR resistance gene
	Multiple cloning site	Start: 5983 End: 6650	pBS replication origin (<i>E. coli</i>)