

Medicago nodules NubG-x cDNA library P02236

Source data

Organism	Medicago nodules
RNA	total RNA

Construction data

Library vector	pPR3-N
Cloned	directional / Sfi I
1st strand synthesis	oligo dT
5' adapter	DSM4

Quality control data

Complexity	1.0 x 10E6 independent clones
Average insert size	1.4 kb
Size range	0.6 - 10.0 kb
% vectors with insert	95%
% inserts > 250 bp	100%

DSM4 adapter	5' AAGCAGTGGTATCAACGCAGAGTGGCCATTACGGCCGGG 3'
oligo dT primer	5' ATTCTAGAGCCGAGGCGGCCGACATGTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTTN 3'

Vector sequences and maps can be found in the support section of <http://www.dualsystems.com>.

Notice to purchaser:

This product is for research use only. This product, or any of its components, may not be transferred for consideration or sold to any third party without the prior written